

The CBL-B Inhibitor, NX-0255, Enhances Human Drug-Enhanced Tumor-Infiltrating Lymphocyte (DeTIL) Expansion and T Cell Function in Full-Scale Runs


Xiaoyan Liang, Xianzhu Wu, Jeevitha Jeevan, Samuel C. Butler, Pranav Murthy, Arthur T. Sands, Michael T. Lotze

Nurix Therapeutics, San Francisco, CA and Pittsburgh, PA, USA

Introduction


- Tumor-infiltrating lymphocyte (TIL) therapy is a promising adoptive cell therapy (ACT) that has demonstrated durable complete responses in patients with metastatic melanoma and several epithelial tumors.^{1,2}
- Obtaining TIL with sufficient quality (increased memory phenotype and reduced exhaustion) during the manufacturing process is a major challenge.^{3,4}
- The E3 ubiquitin ligase, Casitas B-lineage lymphoma B (CBL-B), limits T cell activation by inhibiting T cell receptor signaling in the absence of CD28 (signal 2) engagement^{5,6}. CBL-B also regulates NK, DC, and B cell activation and function.
- Addition of ex vivo NX-0255, a highly potent CBL-B inhibitor (CBL-Bi), significantly increases TIL growth with a favorable phenotype at research scale expansion (Fig. 1).

Figure 1. Overview of Drug-Enhanced TIL Therapy (DeTIL-0255)



Materials and Methods

- We assessed the impact of ex vivo CBL-Bi with NX-0255 on the expansion and phenotype of drug-enhanced TIL (DeTIL), in multiple full-scale processes.
- Various tumor types were assessed, including endometrial (n=2), cervical (n=1), lung (n=1), colon (n=1), and melanoma (n=1) tumors.
- Six full-scale studies were performed in parallel with TIL expanded either solely in the presence of IL-2 (TIL arm), or in the presence of IL-2 and NX-0255 (DeTIL-0255 arm).
- TIL and DeTIL-0255 were harvested following 22 days of expansion and assessed for total cell number, viability, phenotype, and function.



Statistical analysis

All statistical analyses were performed using a two-tailed Wilcoxon signed-rank test. *p<0.05, **p<0.01. Red color represents the DeTIL-0255 arm, and blue color represents the TIL arm.

References

1. Rosenberg SA, et al. Clin Cancer Res. 2011;17:4550–7.
2. Zacharakis N, et al. Nat Med 2018; 24, 724–730.
3. Wang S, et al. BMC Medicine 2021;19:140–7.
4. Krishna S, et al. Science 2020;370:1328–34.
5. Lutz-Micoladoni C, et al. Frontiers in Oncology 2015;5:1–14.
6. Zhou X, et al. Signal Transduction And Targeted Therapy 2021;6:16.

Results

Compared with the TIL arm, the addition of NX-0255 in DeTIL-0255:

- Increases the number of total viable cells in five of six full-scale experiments (Figure 2).
- Significantly increases the total number of CD8⁺ T cells (Figure 3).
- Significantly increases the total number of CD8⁺ & CD4⁺ central memory T cells (Figure 4).
- Significantly increases the number of CD8⁺ CD39⁺ CD69⁺ 'stem-like' memory cells (Figure 5).
- Significantly increases the intracellular expression of IFN- γ in response to PMA and Ionomycin stimulation and granzyme B, in response to α CD3- α CD28 stimulation (Figure 6).

Figure 2. DeTIL-0255 demonstrate increased number of total viable cells in 5 of 6 full-scale manufacturing runs




Figure 3. Addition of NX-0255 increases the total number of CD8⁺ T cells




Figure 4. DeTIL-0255 demonstrate significantly increased CD4⁺ and CD8⁺ central memory T cells compared to TIL




Figure 5. DeTIL-0255 significantly increased CD39-CD69⁺ stem-like CD8⁺ T cell when compared with TIL





Figure 6. DeTIL-0255 demonstrate significant increases in the proportion of CD8⁺ T cells expressing intracellular IFN- γ and granzyme B following stimulation



Conclusions

Manufacturing of DeTIL-0255 is feasible and, in comparison to conventional TIL, demonstrates a superior phenotype that is suitable for clinical evaluation. Accordingly, we have initiated a clinical trial with DeTIL-0255 in patients with gynecologic malignancies (ovarian cancer, endometrial cancer, and cervical cancer), having treated our first patient (NCT05107739).

